

Product Datasheet

MAP4K4 (phospho-Ser629) antibody (orb34671)

Description	Rabbit polyclonal to MAP4K4 (phospho-Ser629)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MAP4K4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	142101
Uniprot ID	O95819
NCBI	NP_001229488.1 , NP_004825.3 , NP_663720.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

