

Product Datasheet

**F(ab')₂ Mouse IgG (H&L) antibody (FITC)
(orb346785)**

Description

F(ab')₂ Mouse IgG (H&L) antibody (FITC)

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

FITC

Tested

FC, FLISA, IF

Applications
Immunogen

Mouse IgG whole molecule

Preservatives

0.01% (w/v) Thimerosal

Form/Appearance

Lyophilized

Concentration

10.0 mg/mL

Storage

Store vial at 4° C prior to restoration.
For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.
Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

Isotype

IgG F(ab')₂
Clonality

Polyclonal

Purity

This product is a F(ab')₂ fragment of an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution Range

FLISA: 1:10,000 - 1:50,000, FC: 1:500 -