

Product Datasheet

PIK3CG (phospho-Ser1100) antibody (orb34707)

Description Rabbit polyclonal to PIK3CG (phospho-Ser1100)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target PIK3CG

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 126454

Uniprot ID **P48736**

NCBI **NP_002640.2, NP_001269356.1, NP_001269355.1**

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

