

Product Datasheet

LRRC15 antibody (HRP) (orb34710)

Description	LRRC15 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 15 protein (301-538AA)
Target	LRRC15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TF66
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.