

Product Datasheet

AKT3 (phospho-Ser472) antibody (orb34711)

Description	Rabbit polyclonal antibody to AKT3 (phospho-Ser472)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	AKT3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55775
Uniprot ID	Q9Y243
NCBI	NP_005456.1 , NP_859029.1 , NP_001193658.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

