

Product Datasheet

Golden Syrian Hamster IgG (H&L) antibody (Alkaline Phosphatase) (orb347125)

Description

Golden Syrian Hamster IgG (H&L)

Species/Host

Rabbit

Reactivity

Golden Hamster

Conjugation

AP

Tested

ELISA, IHC, WB

Applications
Immunogen

Golden Syrian Hamster IgG, whole molecule

Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Liquid (sterile filtered)

Concentration

1.0 mg/mL

Storage

Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.

Note

For research use only

Application notes

Anti-Golden Syrian Hamster IgG Alk Phos conjugate is suitable for immunoblotting (western or dot blot), ELISA, immunoelectron microscopy and immunohistochemistry as well as other antibody-based enzymatic assays requiring lot-to-lot consistency.

Isotype

IgG

Clonality

Polyclonal

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Golden Syrian Hamster IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Rabbit Serum, Golden Syrian Hamster IgG, and