

Product Datasheet

TGFBRAP1 antibody (HRP) (orb34720)

Description	TGFBRAP1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Transforming growth factor-beta receptor-associated protein 1 protein (401-600AA)
Target	TGFBRAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WUH2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.