

Product Datasheet

PIK3CD (phospho-Tyr485) antibody (orb34735)

Description Rabbit polyclonal to PIK3CD (phospho-Tyr485)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT

Immunogen Synthetic Peptide

Target PIK3CD

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 119479

Uniprot ID [O00329](#)

NCBI [NP_005017.3](#)

Dilution Range DOT: 1:500

Expiration Date 12 months from date of receipt.

