

Product Datasheet

ST8SIA4 antibody (orb34747)

Description	ST8SIA4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human CMP-N-acetylneuraminate-poly-alpha-2,8-sialyltransferase protein (21-168AA)
Target	ST8SIA4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92187
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.