

Product Datasheet

ND-75 antibody (Biotin) (orb34770)

Description	ND-75 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Drosophila melanogaster NADH-ubiquinone oxidoreductase 75 kDa subunit, mitochondrial protein (28-315AA)
Target	ND-75
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q94511
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.