

## Product Datasheet

**Rat IgG IgA IgM (H&L) antibody (Peroxidase)  
(orb347791)**

**Description**

Rat IgG IgA IgM (H&amp;L) antibody

**Species/Host**

Goat

**Reactivity**

Rat

**Conjugation**

HRP

**Tested**

ELISA, IHC, WB

**Applications**
**Immunogen**

Rat IgG, IgA and IgM whole molecules

**Preservatives**

0.01% (w/v) Thimerosal

**Form/Appearance**

Lyophilized

**Concentration**

2.0 mg/mL

**Storage**

Store vial at 4° C prior to restoration.  
For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.  
Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

**Note**

For research use only

**Application notes**

Anti-Rat IgG IgA IgM Antibody Peroxidase Conjugated has been tested by ELISA and is suitable for immunoassays where reactivity to different classes of rat immunoglobulins is desired. Antibody has been tested by ELISA, western blot, and immunohistochemistry. Optimal concentrations in immunoassays should be determined by the researcher.

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Anti-Rat IgG IgA IgM (H&L) Antibody, Peroxidase Conjugated, was prepared from polyspecific antiserum by immunoaffinity chromatography using Rat antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive&lt;br&gt;Research Triangle Park&lt;br&gt;Durham, North Carolina&lt;br&gt;27709, United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558