



Product Datasheet

DAZAP1 antibody (orb34783)

Description	DAZAP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human DAZ-associated protein 1 protein (1-300AA)
Target	DAZAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96EP5
Dilution Range	WB:1:2000-10000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), J...



Immunofluorescence analysis of HeLa cell...