



Product Datasheet

ZDHHC12 antibody (orb34797)


Description	ZDHHC12 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC12 protein (121-197AA)
Target	ZDHHC12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96GR4
Dilution Range	WB:1:500-2000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of human th...