

Product Datasheet

Donkey IgG (H&L) antibody (Alkaline Phosphatase) (orb347995)

Description

Donkey IgG (H&L) antibody (Alkaline

Species/Host	Goat
Reactivity	Donkey
Conjugation	AP
Tested Applications	DOT, ELISA, IHC, WB
Immunogen	Donkey IgG whole molecule
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Liquid (sterile filtered)
Concentration	1.0 mg/mL
Storage	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.
Note	For research use only
Application notes	Anti-Donkey IgG Alkaline Phosphatase Antibody has been tested by ELISA and dot blot and is suitable for immunoblotting (western or dot blot), ELISA and immunohistochemistry as well as other phosphatase-antibody based enzymatic assays requiring lot-to-lot consistency.
Isotype	IgG
Clonality	Polyclonal
Purity	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Donkey IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Goat Serum, Donkey IgG, and Donkey Serum.
Dilution Range	ELISA: 1:10,000 - 1:40,000, IHC: 1:200 - 1:1,000, WB: 1:500 - 1:2,500