



Product Datasheet

CAPNS2 antibody (orb34834)

Description	CAPNS2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Calpain small subunit 2 protein (1-248AA)
Target	CAPNS2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96L46
Dilution Range	WB:1:500-2000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart tis...



Immunofluorescence analysis of A549 cell...