

## Product Datasheet

### Fab Guinea Pig IgG (H&L) antibody (Biotin) (orb348437)

**Description**

Fab Guinea Pig IgG (H&L) antibody

|                            |                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Goat                                                                                                                                                                                                                                                                                                                                                  |
| <b>Reactivity</b>          | Guinea pig                                                                                                                                                                                                                                                                                                                                            |
| <b>Conjugation</b>         | Biotin                                                                                                                                                                                                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                        |
| <b>Immunogen</b>           | Guinea Pig IgG whole molecule                                                                                                                                                                                                                                                                                                                         |
| <b>Preservatives</b>       | 0.01% (w/v) Sodium Azide                                                                                                                                                                                                                                                                                                                              |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                                                                                                                                           |
| <b>Concentration</b>       | 1.0 mg/mL                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>             | Store vial at 4° C prior to restoration.<br>For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.<br>Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                 |
| <b>Application notes</b>   | Suitable for immunoblotting, ELISA, immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.                                                            |
| <b>Isotype</b>             | IgG Fab                                                                                                                                                                                                                                                                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>              | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Guinea Pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a                                         |

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive<br>Research Triangle Park<br>Durham, North Carolina<br>27709. United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558