



Product Datasheet

PSKH2 antibody (orb34848)

Description	PSKH2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Serine/threonine-protein kinase H2 protein (63-320AA)
Target	PSKH2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96QS6
Dilution Range	WB:1:500-1000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cells(lane...)



Immunohistochemical staining of human gastric tissue (ga...)