

Product Datasheet

Alkaline Phosphatase (orb348523)

Description	Alkaline Phosphatase
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None. 6mM Tris Chloride, 6mM Magnesium Chloride, 0.12mM Zinc Chloride; pH 7.6
Form/Appearance	Liquid (sterile filtered)
Concentration	18.2 mg/ml
Storage	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity.
Note	For research use only
Purity	This product was prepared from chromatographically pure alkaline phosphatase (calf intestinal mucosa).
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.