

Product Datasheet

1.0 M Tris Hcl Ph 8.0 (orb348529)

Description

1.0 M Tris HCl Ph 8.0

Conjugation

Unconjugated

Preservatives

Preservative: None. Stabilizer: None. See application note.

Form/Appearance

Liquid (sterile filtered)

Concentration

1.0 M

Storage

Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.

Note

For research use only

Application notes

Tris solution is suitable for molecular immunology, molecular biology and nucleic acid methodologies, when applicable. Visit our newly expanded web site at www.rockland.com for methods using this and other buffers. This product is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. This buffer consists of 1.0 M Tris HCl at pH 8.0. Meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade (WFI) water.

Purity

Tris HCl 1.0 M at pH 8.0 is meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade (WFI) water.

Hazard Information

Non-Toxic

Expiration Date

12 months from date of receipt.