

Product Datasheet

Ultrapure Sterile Water (orb348534)

Description	Ultrapure Sterile Water
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less.
Note	For research use only
Application notes	Ultra Sterile Water is suitable for molecular immunology, molecular biology and nucleic acid methodologies, when applicable. Visit our web site at www.rockland.com for methods using this and other buffers.
Purity	UltraPure Sterile Water consists of highly polished pharmaceutical grade deionized water processed using sterile conditions This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.