

Product Datasheet

Ultrapure Sterile Water (orb348535)

| | |
|---------------------------|--|
| Description | Ultrapure Sterile Water |
| Conjugation | Unconjugated |
| Preservatives | Preservative: None. Stabilizer: None |
| Form/Appearance | Liquid (sterile filtered) |
| Concentration | 1X |
| Storage | Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. |
| Note | For research use only |
| Application notes | Suitable for molecular immunology, molecular biology and nucleic acid methodologies, when applicable. Visit our web site at www.rockland.com for methods using this and other buffers. |
| Purity | UltraPure Sterile Water consists of highly polished pharmaceutical grade deionized water processed using sterile conditions This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria. |
| Hazard Information | Non-Toxic |
| Expiration Date | 12 months from date of receipt. |