

Product Datasheet

Molecular Biology Grade Ultrapure Water (orb348537)

Description

Molecular Biology Grade Ultrapure

Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less.
Note	For research use only
Application notes	Molecular Biology Grade Ultrapure water is suitable for molecular immunology, molecular biology and nucleic acid methodologies, when applicable. Visit our web site at www.rockland.com for methods using this and other buffers.
Purity	Molecular Biology Grade UltraPure Water consists of highly polished pharmaceutical grade deionized water treated with diethyl pyrocarbonate (DEPC). This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.