

Product Datasheet

HIPK2 (phospho-Tyr361) antibody (orb34854)

Description	Rabbit polyclonal antibody to HIPK2 (phospho-Tyr361) belongs to the protein kinase super-family. CMG...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	HIPK2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	130966
Uniprot ID	Q9H2X6
NCBI	NP_001106710.1 , NP_073577.3
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

