

Product Datasheet

10% Sodium Dodecyl Sulfate Sds (orb348542)

Description	10% Sodium Dodecyl Sulfate Sds
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None. See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	10% (w/v)
Storage	Store container at room temperature (18° to 26° C) prior to opening.
Note	For research use only
Application notes	This product is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. Product consists of 10% (w/v) Sodium Dodecyl Sulfate (SDS). Meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade deionized water.
Purity	10% SDS was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.