

## Product Datasheet

### Get Buffer (orb348563)

<b>Description</b>	Get Buffer
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	Preservative: None. Stabilizer: None. See application note.
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Storage</b>	Store container at 4° C prior to opening. Use aseptic technique to maintain sterility when opening product.
<b>Note</b>	For research use only
<b>Application notes</b>	This product is a 1X concentrated stock solution. No dilution is required. This buffer consists of 50 mM Glucose, 10 mM EDTA, and 25 mM TRIS Cl at a pH of 8.0. Meticulously prepared using ultra pure reagents dissolved in pharmaceutical grade deionized water.
<b>Purity</b>	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
<b>Hazard Information</b>	Non-Toxic
<b>Expiration Date</b>	12 months from date of receipt.