

Product Datasheet

1X Tes Ph 7.8 (orb348575)

Description	1X Tes Ph 7.8
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None. See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
Note	For research use only
Application notes	This product is ready-to-use as a working 1X solution and requires no further dilution. This buffer consists of 10 mM Tris HCl, 100 mM NaCl, and 1.0 mM EDTA at a pH of 7.8. Meticulously prepared using ultra pure reagents dissolved in DEPC treated pharmaceutical grade deionized water.
Purity	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.