

## Product Datasheet

### Pbs Fish Gel Concentrate 10X (orb348587)

**Description**

Pbs Fish Gel Concentrate 10X

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, WB

**Preservatives**

See application note.

**Form/Appearance**

Liquid (sterile filtered)

**Concentration**

10X

**Storage**

Store container at 4° C before opening.  
Protect from moisture and light. No special shipping conditions or precautions are required.

**Note**

For research use only

**Application notes**

This product is a 10X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 10X concentrate with 9 parts distilled-deionized water or equivalent. 10X PBS Fish Gel Concentrate consists of 0.2M potassium phosphate, 1.5M sodium chloride and Fish Gelatin at pH 7.2. A proprietary combination of stabilizers and preservatives are used that are azide and mercury free.

**Purity**

This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.

**Hazard Information**

Non-Toxic

**Dilution Range**

ELISA: User Optimized, WB: User Optimized

**Expiration Date**

12 months from date of receipt.