

Product Datasheet

Tbs Fish Gel Concentrate 10X (orb348589)

Description

Tbs Fish Gel Concentrate 10X

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Preservatives

See application note.

Form/Appearance

Liquid (sterile filtered)

Concentration

10X

Storage

Store container at 4° C before opening.
Protect from moisture and light. No special shipping conditions or precautions are required.

Note

For research use only

Application notes

This product is a 10X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 10X concentrate with 9 parts distilled-deionized water or equivalent. 10X TBS Fish Gel Concentrate consists of 1.0 M Tris hydrochloride, 1.5 M sodium chloride and Fish Gelatin at pH 7.5. A proprietary combination of stabilizers and preservatives are used that are azide and mercury free.

Purity

This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.

Hazard Information

Non-Toxic

Dilution Range

ELISA: User Optimized, WB: User Optimized

Expiration Date

12 months from date of receipt.