

Product Datasheet

Tbs Fish Gel 10X (orb348590)

Description	Tbs Fish Gel 10X
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Preservatives	See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	10X
Storage	Store container at 4° C before opening. Protect from moisture and light. No special shipping conditions or precautions are required.
Note	For research use only
Application notes	This product is a 10X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 10X concentrate with 9 parts distilled-deionized water or equivalent. 10X TBS Fish Gel Concentrate consists of 1.0 M Tris hydrochloride, 1.5 M sodium chloride and Fish Gelatin at pH 7.5. A proprietary combination of stabilizers and preservatives are used that are azide and mercury free.
Purity	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Dilution Range	ELISA: User Optimized, WB: User Optimized
Expiration Date	12 months from date of receipt.