

Product Datasheet

Ponceau S Solution (orb348592)

Description	Ponceau S Solution
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid
Concentration	1X
Storage	Store container at room temperature (20° to 26° C).
Note	For research use only
Application notes	Store at room temperature. The product is at working concentration, no further dilution is necessary. The limit of detection for this stain is 200 nanograms of protein after separation by electrophoresis in polyacrylamide gels and transferred to nitrocellulose or PVDF membranes. Ponceau S is not suitable for use with nylon membranes.
Purity	Ponceau S Staining Solution contains Ponceau S in a specially formulated acetate buffer for optimal staining results. Ponceau S is a ready-to-use stain that can be used to quickly and reversibly stain proteins on blotting membranes after Western blot transfer. Staining of Nylon membranes is not recommended.
Hazard Information	Non-Toxic
Dilution Range	WB: ready to use
Expiration Date	12 months from date of receipt.