

Product Datasheet

10X Tris-Borate-Edta Tbe Buffer (orb348612)

Description	10X Tris-Borate-Edta Tbe Buffer
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None. See application note.
Form/Appearance	Liquid (sterile filtered)
Concentration	10X
Storage	Store container at room temperature (18° to 26° C) prior to opening.
Note	For research use only
Application notes	This product is a concentrated stock solution and should be diluted appropriately with distilled, deionized water or equivalent to its final working concentration. 10X Tris-Borate EDTA (TBE) consists of 0.9 M Tris-Borate, 0.01 M EDTA at a pH of 8.3. Meticulously prepared using ultra pure reagents dissolved in highly polished pharmaceutical grade deionized water.
Purity	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.