

## Product Datasheet

### Normal Saline (orb348615)

<b>Description</b>	Normal Saline
<b>Conjugation</b>	Unconjugated
<b>Preservatives</b>	Preservative: None. Stabilizer: None. 0.15 M Sodium Chloride
<b>Form/Appearance</b>	Liquid (sterile filtered)
<b>Concentration</b>	1X
<b>Storage</b>	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
<b>Note</b>	For research use only
<b>Application notes</b>	Saline solutions are available in various formulations for different purposes. Salines are also used in cell biology, molecular biology, and biochemistry experiments where ever normal salines are needed.
<b>Purity</b>	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
<b>Hazard Information</b>	Non-Toxic
<b>Expiration Date</b>	12 months from date of receipt.