

Product Datasheet

Borate Buffered Saline (Bbs) (orb348616)

Description	Borate Buffered Saline (Bbs)
Conjugation	Unconjugated
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store container at room temperature (18° to 26° C) prior to opening. If desired, the solution may be stored at 4° C or less. Some salts may precipitate out of solution at lower temperature. Allow buffer to equilibrate to room temperature (18° to 26° C) to restore solubility of some salts.
Note	For research use only
Application notes	BBS is isotonic and has a strong bactericidal effect. It can be used to dilute solutions for coating procedures. Additives such as Tween 20 and Blotto can be used to add functionality as a washing buffer or blocking buffer.
Purity	This product was aseptically filtered through a Millipore 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.