

Product Datasheet

Blocking Buffer (2X) For Fluorescent Western Blotting (orb348639)

Description

Blocking Buffer (2X) For Fluorescent

Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Thimerosal is added as an antimicrobial agent.
Form/Appearance	Liquid (sterile filtered)
Concentration	2X
Storage	Store blocking buffer at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	This product is a 2X concentrated stock solution. Prepare a 1X working solution by diluting 1 part 2X concentrate with 1 parts TBS or equivalent. This buffer was meticulously prepared using ultra pure reagents dissolved in pharmaceutical grade water (WFI) and consists of a proprietary protein formulation in TRIS buffered saline at pH 7.6 with thimerosal added as an antimicrobial agent.
Purity	Blocking buffer is specifically formulated to achieve superior reproducible Western blotting images using this system.
Hazard Information	Non-Toxic
Dilution Range	WB: User Defined
Expiration Date	6 months from date of receipt.