

Product Datasheet

Blocking Buffer For Fluorescent Western Blotting Thimerisol (orb348641)

Description

Blocking Buffer For Fluorescent

Conjugation	Unconjugated
Tested Applications	WB
Preservatives	0.01% (w/v) Gentamicin Sulfate
Form/Appearance	Liquid (sterile filtered)
Concentration	1X
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE.
Note	For research use only
Application notes	Optimum performance is achieved using this product undiluted. However, this blocking solution can often be diluted 1:1 in TBS without a significant loss of performance. This buffer was prepared using ultra pure reagents dissolved in pharmaceutical grade water (WFI) and consists of a proprietary protein formulation in TRIS buffered saline at pH 7.6.
Purity	Blocking buffer was prepared using ultra pure reagents dissolved in pharmaceutical grade water (WFI) and consists of a proprietary protein formulation in TRIS buffered saline at pH 7.6.
Hazard Information	Non-Toxic
Dilution Range	WB: User Defined
Expiration Date	6 months from date of receipt.