




Product Datasheet

Revitablot Western Blot Stripping Buffer (orb348648)

Description	Revitablot Western Blot Stripping Buffer
Conjugation	Unconjugated
Tested Applications	WB
Preservatives	Preservative: None. Stabilizer: None
Form/Appearance	Liquid
Concentration	1X
Storage	Store container at room temperature (18° to 26° C) prior to opening.
Note	For research use only
Application notes	Revitablot Western Blot Stripping Buffer is ready-to-use as a working 1X solution and requires no further dilution.
Purity	Revitablot Western Blot Stripping Buffer contains solutions in a proprietary combination to enhance the removal of bound antibodies from western blot membranes for repeated use. The proprietary formulation of the solution ensures high stripping efficiency with low backgrounds.
Hazard Information	Non-Toxic
Dilution Range	WB: Use undiluted.
Expiration Date	12 months from date of receipt.



Western
Blot
Stripping
Buffer.
Lane 1: 7...



Western
Blot
Stripping
Buffer
using
Fluo...