

Product Datasheet

Tmb Elisa Peroxidase Substrate (orb348653)

Description	Tmb Elisa Peroxidase Substrate
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Proprietary buffer with enhancer
Form/Appearance	Liquid - clear to very light blue colored liquid
Concentration	1X
Storage	Store container at 4° C prior to opening. Protect from moisture and light. No special shipping conditions or precautions are required.
Note	For research use only
Application notes	TMB ELISA Peroxidase Substrate will produce a soluble blue end product read at 370 nm or 655 nm. TMB ELISA Peroxidase Substrate incubation time will vary depending on the assay conditions. No dilutions are required, TMB ELISA Peroxidase Substrate comes ready to use.
Purity	pH: 4.0 +/- 0.05 Stability @ 18° to 26°C: PASS Stability @ 4°C: PASS QC Raw Material: PASS Absorbance check of final product: PASS Performance Data per ELISA: PASS
Hazard Information	Non-Toxic
Dilution Range	ELISA: 1X
Expiration Date	Please enquire.