

Product Datasheet

Tmb Prestained Red Elisa Peroxidase Substrate (orb348655)

Description	Tmb Prestained Red Elisa Peroxidase Substrate
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Proprietary buffer with enhancer
Form/Appearance	Liquid - light red colored
Concentration	1x
Storage	Store container at 4° C prior to opening. Protect from moisture and light. No special shipping conditions or precautions are required.
Note	For research use only
Application notes	TMB Prestained Red ELISA Peroxidase Substrate produces a soluble blue end product of ELISA, read at 370nm or 655nm.
Purity	pH: 4.0 +/- 0.2 Stability @ 18° to 26°C: PASS Stability @ 4°C: PASS QC Raw Material: PASS Absorbance check of final product: PASS Performance Data per ELISA: PASS
Hazard Information	Non-Toxic
Dilution Range	ELISA: 1x
Expiration Date	Please enquire.