

Product Datasheet

Tmb Membrane Peroxidase Substrate (orb348656)

Description	Tmb Membrane Peroxidase Substrate
Conjugation	Unconjugated
Tested Applications	DOT, IHC, WB
Preservatives	Proprietary
Form/Appearance	Liquid - clear
Concentration	1X
Storage	Store container at 4° C prior to opening. Protect from moisture and light. No special shipping conditions or precautions are required.
Note	For research use only
Application notes	TMB Membrane Peroxidase Substrate comes ready to use. No dilutions are required. TMB Membrane Peroxidase Substrate will produce an insoluble blue-green end product. Positive signal deposits blue-green color on membrane or surface of sample material. TMB Membrane Peroxidase Substrate incubation time will vary depending on the assay conditions.
Hazard Information	Non-Toxic
Dilution Range	IHC: 1X, WB: 1X
Expiration Date	12 months from date of receipt.