

Product Datasheet

Protein A (orb348747)

Description	Protein A
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	Lyophilized
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Protein A is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Protein A is suitable for use as an antigen, as a control or standard in assays, for conjugation and for most other immunological methods. Protein A contains only IgG binding domains, ensuring maximum specific IgG binding.
Purity	Protein A is chromatographically pure and shows predominantly a single band by SDS-PAGE. Greater than 95% of the A280 material binds to Human IgG.
Uniprot ID	P38507
Hazard Information	Non-Toxic
Dilution Range	ELISA: User Optimized, WB: User Optimized
Expiration Date	6 months from date of receipt.