

Product Datasheet

Protein A Fluorescein Conjugated (orb348748)

Description

Protein A Fluorescein Conjugated

Conjugation

FITC

Tested

DOT, FC, FLISA, IF, WB

Applications
Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Lyophilized

Concentration

1.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Protein A Fluorescein Conjugated has been tested by dot blot and western blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.

Purity

This product was prepared from chromatographically pure Protein A. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein and anti-Protein A.

Uniprot ID
P38507
Hazard

Non-Toxic

Information
Dilution Range

FLISA: 1:10,000 - 1:50,000, FC: 1:500 - 1:2,500, IF: 1:1,000 - 1:5000

Expiration Date

6 months from date of receipt.