

Product Datasheet

C4orf33 antibody (FITC) (orb34875)

Description	C4orf33 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human UPF0462 protein C4orf33 protein (1-199AA)
Target	C4orf33
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N1A6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.