

Product Datasheet

Protein A Biotin Conjugated (orb348752)

Description	Protein A Biotin Conjugated
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Lyophilized
Concentration	1.0 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Protein A Biotin Conjugated is suitable as a detection agent for primary antibodies that are of the IgG isotype.
Purity	This product was prepared from chromatographically pure Protein A. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin and anti-Protein A.
Uniprot ID	P38507
Hazard Information	Non-Toxic
Dilution Range	ELISA: 1:20,000 - 1:200,000, IHC: 1:1,000 - 1:5,000, WB: 1:10,000 - 1:40,000
Expiration Date	6 months from date of receipt.