

Product Datasheet

Sepharose Protein G (orb348763)

Description	Sepharose Protein G
Conjugation	Sepharose
Tested Applications	IP
Form/Appearance	Suspension of agarose beads
Storage	Store SEPHAROSE PROTEIN G at 4° C prior to opening. DO NOT FREEZE. Supplied in 20% Ethanol/Water
Note	For research use only
Application notes	Protein G resin may also be used to pull down antibody:antigen complexes in immunoprecipitation experiments. This product is suitable for use at 20 µL per immunoprecipitation reaction.
Purity	Sepharose Protein G binds to all IgG subclasses from human, monkey, mouse and rat species. It also binds to total IgG from guinea pig, rabbit, goat, cow, sheep, and horse.
Uniprot ID	P19909
Hazard Information	Non-Toxic
Dilution Range	IP: User Optimized
Expiration Date	6 months from date of receipt.