

Product Datasheet

**Streptavidin Atto 532 Conjugated [Out of stock]
(orb348775)**

Description

Streptavidin Atto 532 Conjugated

Conjugation

ATTO-532

Tested

DOT, FLISA, IF, WB

Applications
Preservatives

0.01% (w/v) Sodium Azide

Form/Appearance

Lyophilized

Concentration

1.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

The emission spectra for this ATTO conjugate matches the principle output wavelengths of most common fluorescence instrumentation.

Purity

STREPTAVIDIN ATTO 532 was prepared from chromatographically purified Streptavidin. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Streptavidin.

Hazard
Information

Non-Toxic

Dilution Range

FLISA: >1:20,000, IF: >1:5,000, WB: >1:10,000

Expiration Date

12 months from date of receipt.