


Product Datasheet

Mouse TIGIT Protein (orb348797)

Description	T-cell immunoreceptor with Ig and ITIM domains (TIGIT) is also known as V-set and immunoglobulin dom...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIGIT
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 14.8 kDa. The protein migrates as 18-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001139797.1
Purity	95%
Source	Mouse TIGIT, His Tag (orb348797) is expressed from human 293 cells (HEK293). It contains AA Gly 26 - Thr 143 (Accession # NP_001139797.1).
MW	14.8 kDa
NCBI	NP_001139797.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse TIGIT
protein