



Product Datasheet

Human ROR2 / NTRKR2 Protein (orb348799)

Description	Human ROR2, His Tag is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ROR2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 43.2 kDa. The protein migrates as 55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_004551.1
Purity	95%
Source	Human ROR2, His Tag (orb348799) is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Accession # NP_004551.1).
MW	43.2 kDa
Uniprot ID	A1L4F5
NCBI	NP_004551.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human ROR2 protein



Line graph illustrates about the Ag-Ab r...