

Product Datasheet

Human HVEM / TNFRSF14 Protein (orb348800)

Description

Human HVEM, His Tag (HVM-H52E9) is expressed

kDa	M	R
116.0		

39 - ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

HVEM

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.2 kDa. The protein migrates as 30-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_003811.2

Purity

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human HVEM, His Tag (orb348800) is expressed from human 293 cells (HEK293). It contains AA Leu 39 - Val 202 (Accession # Q92956-1).

MW

19.2 kDa

Uniprot ID


[Q92956](#)

NCBI

[NP_003811.2](#)

Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human HVEM protein