


Product Datasheet

Mouse Transferrin R2 Protein (orb348801)

Description	Mouse TFR2, Mouse Fc Tag is expressed from human 293 cells (HEK293). It contains AA Arg 103 - Phe 79...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Transferrin R2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	N-mFc(mIgG1)
Note	For research use only
Application notes	This protein carries a mouse IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 101.1 kDa. The protein migrates as 116 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_056614.3
Purity	95%
Source	Mouse TFR2, Mouse Fc Tag (orb348801) is expressed from human 293 cells (HEK293). It contains AA Arg 103 - Phe 798 (Accession # Q9JKX3-1).
MW	101.1 kDa
Uniprot ID	Q9JKX3
NCBI	NP_056614.3
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse
Transferrin
R...