

Product Datasheet

Human Nectin-2 / CD112 Protein (orb348803)

Description Human Nectin-2, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 32 - Leu 360 ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Nectin-2

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.4 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_002847.1

Purity 95%


Source Human Nectin-2, His Tag (orb348803) is expressed from human 293 cells (HEK293). It contains AA Gln 32 - Leu 360 (Accession # Q92692-2).

MW 37.4 kDa

Uniprot ID [Q92692](#)

NCBI [NP_002847.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Nectin-2
prot...